

ASPIRATION ACCESSORIES

CANISTER COLLECTION SYSTEMS

When size matters, we offer three great choices.

BAXTER™ 1500CC SYSTEM

Baxter™ Canister
Reusable 1500cc 20-5153-00

Baxter™ Liner
Disposable 1500cc 20-5154-00 (Case of 10)

SORENSEN™ 1900CC SYSTEM

Sorenson™ Canister
Reusable 1900cc 20-5157-00

Sorenson™ Liner
Disposable 1900cc 20-5158-00 (Case of 10)

QUICK-FIT™ 3000CC SYSTEM

Quick-Fit™ Canister
Reusable 3000cc 20-5156-00

Quick-Fit™ Liner
Disposable 3000cc 20-5155-00 (Case of 10)

BAXTER™ CANISTER HOLDER

Baxter Canister Holder fits on top panel of the aspirator to hold 1500cc Baxter™ canisters in place.
Part# 20-5159-01

SORENSEN™ CANISTER HOLDER

Sorenson Canister Holder fits on top panel of the aspirator to hold 1900cc Sorenson™ canisters in place.
Part# 20-5159-00

QUICK-FIT™ CANISTER HOLDER

Quick-Fit Canister Holders fits on top panel of the aspirator to hold 3000cc Quick-Fit™ canisters in place.
Part# 20-5159-03



TUBING CHOICES

Wells Johnson offers the most complete selection of infiltration products in the world!
Choose single spike, "Y" connector, or large bore. Ready to use - no assembly required

ASPIRATION TUBING

STANDARD ULTRA-LIMP™

Part# 24-5103-00 (Case of 10)

9ft long for use with Baxter canister

Will withstand 29.5" vacuum (will not collapse)



MB ULTRA-LIMP™

Part# 24-5104-00 (Case of 10)

Length 9ft

Will withstand 29.5" vacuum (will not collapse)

Also fits QuickFit system



INFILTRATION TUBING

HVP REINJECTION

Part# 24-6037-00 (case of 10)

Features:

Proprietary tubing

Designed to only be used with the HVP™ unit

144" long



Y INFILTRATION


Part# 24-6009-00 (case of 10)

Features:

Silicone sleeve for use in infiltration

170" total length

Dual spikes for non-stop infiltration



LARGE BORE INFILTRATION

Part# 24-6028-00 (case of 10)

Features:

Single spike

Works with infiltration pump

188" total length

Large bore diameter for increased flow



CB INFILTRATION


Part# 24-6008-00 (case of 10)

Features:

Single spike

Works with infiltration pump

156" total length



POWER ASSISTED TUBING ADAPTER

Power-assisted tubing adapts handpiece to standard Wells Johnson aspiration tubing.

24-5109-44

Reusable. Sold individually.



ACCESSORY PRODUCTS FILTERS AND ACCESSORIES OTHER ASPIRATOR SUPPLIES

BACTERIAL FILTERS

HIGH FLOW

Part# 20-5201-00

.2 micron , maximum use 24 hours

Preventive Maintenance:

Change every 24 hours of use

(Case of 10)



OVERFLOW TRAP ASSEMBLY

ASPIRATOR ACCESSORY

Part# 17-5100-00

Includes jar, lid and overflow ball.



OVERFLOW BALL

Part# 16-5107-00

Replacement for trap assembly.



MUFFLERS

Part# 16-5073-00

Replacement for

Hercules and Whisperator models

Preventive Maintenance:

Change every 80 hours of use.



Muffler replacement for
Wells Johnson Aspirator II
and II+ models.

Part# 16-5109-00

Preventive Maintenance:

Change every 3 months.

